
LEGEND:

- Luer Lock Needle Nut, 2. Needle Nut Seal Ring, 3. Fibre Seal Washer, 4. Cylinder (50mL), 5. Cylinder (30mL), 6. Cylinder (25mL), 7. Delivery Cage Seal Ring, 8. Luer Lock Mount Delivery Cage, 9. Push Rod Guide, 10. Lever, 11. Slide Linkage Pin, 12. Repeater Pawl, 13. Pawl Return Spring, 14. Lever Arm 15. Lever Arm Pivot Pin, 16. Roller, 17. Return Spring, 18. Handle (Left), 19. Needle Magazine, 20. Dose Adjustor (50mL), 21. Dose Adjustor (30mL), 22. Dose Adjustor (25mL), 23. Dose Adjustor Locking Sleeve, 24. Release Locking Ring, 25. Ratchet Plate Spring, 26. Ratchet Plate, 27. Screw (5 Per), 28. Handle (Right), 29. Piston Seal Ring (50mL), 30. Piston Seal Ring (30mL), 31. Piston Seal Ring (25mL), 32. Push Rod Assembly (50mL), 33. Push Rod Assembly (30mL), 34. Push Rod Assembly (25mL), 35. Push Rod End Knob.

SPARE PARTS AVAILABLE:

KIT No.	PART NAME	INCLUDES ITEM No.
WX999	25mL Cylinder and Delivery Cage Assembly	1,2,3,6,7,8,31.
WX1119	30mL Cylinder and Delivery Cage Assembly	1,2,3,5,7,8,30.
WX1111	50mL Cylinder and Delivery Cage Assembly	1,2,3,4,7,8,29.

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Please quote Kit No. and Part Name or instrument code number (PAS802,PAS803 or PAS1078) and item number when ordering spare parts.

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No. de KIT	NOM DE LA PIECE	COMPRED ELEMENT ref.
WX999	Kit soupape / cylindre 25 ml	1,2,3,6,7,8,31.
WX1119	Kit soupape / cylindre 30 ml	1,2,3,5,7,8,30.
WX1111	Kit soupape / cylindre 50 ml	1,2,3,4,7,8,29.

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Semi-Auto Injector

WITH LUER LOCK

njphillips.com
Our expertise is in your hands.

The 25mL, 30mL and 50 mL Semi Automatic Injectors are designed for administration of most vaccines and injectable solutions to livestock.

As components in these instruments may be affected by solvents in some commonly used farm chemicals no responsibility will be accepted by the manufacturer should these instruments be used with such products.

INSTRUCTIONS FOR USE
A. Sterilization

It is essential that this instrument and a supply of needles be thoroughly sterilized before each use. A common method of sterilization is as follows:

- Draw hot water into the cylinder.
- Suspend complete instrument in a container of water and boil together with the needles, for 10 to 20 minutes.
- Remove instrument from container, wrap cloth around the handle and pump dry.

NOTE: Suspending the handpiece not only makes it easier to remove but also prevents damage should the container boil dry. Chemical sterilization with antiseptic solutions is sometimes practised and in such circumstances the recommendations of the chemical manufacturer should be followed. Do not attempt to sterilize by autoclaving.

B. Filling

Before filling from a rubber-sealed glass bottle, ensure that the release locking ring (24) is positioned fully back to disengage the ratchet, allowing the push rod to move freely in either direction. Withdraw push rod fully, expel all air and set injector to the required maximum capacity. Fit needle to syringe and pierce rubber bung on inverted bottle. Introduce air into bottle by pushing the push rod end knob (35). Then, keeping bottle inverted, draw liquid into the cylinder. Repeat this procedure until the syringe is fully charged. Presence of excessive air bubbles indicates an air leak in the instrument or incorrect filling technique.

WARNING: Care must be taken to ensure the liquid does NOT come into contact with any part of the operators body. Chemicals may cause injury to the operator.

C. Dose Setting

Set the dose by turning the dose adjustor (20, 21 or 22) until the required dose is aligned with the arrow on the top of the instrument.

CARE AND MAINTENANCE AFTER USE

- Remove all traces of material from syringe by flushing thoroughly with a hot water detergent mix. Follow with clean, warm water.
- Lubricate piston by unscrewing the cylinder and placing a small quantity of NJ Phillips Lubricant, onto the 'O' ring.
- Place a few drops of NJ Phillips Lubricant on all moving parts.
- When replacing the cylinder, make sure that the push rod is fully back for ease of location, otherwise damage to the cylinder could occur.

DO NOT store your applicator or feedtube full of product. Clean as per the "Care and Maintenance" instructions.

REPLACING PISTON RING

Should it become necessary to replace the piston seal ring, the subsequent procedure must be followed.

- Pull push rod fully back.
- Unscrew cylinder from handle.
- Remove worn 'O' ring from piston, taking care not to damage the groove.
- Lubricate new ring with NJ Phillips Lubricant and fit to piston.
- Replace cylinder as directed in Care and Maintenance.

IMPORTANT
Do not attempt to dismantle the dose adjustor assembly as this should never be necessary.

ABN 36 000 082 002
PAS802 (25mL), PAS803 (50mL), PAS1078 (30mL)

Postal Address: Locked Bag 8
Gosford NSW 2250 AUSTRALIA
Telephone: (02) 4340 2044
Fax: (02) 4340 1991
Email: njp1@njphillips.com.au
www.njphillips.com

FABRIQUE EN AUSTRALIE
Dirección Postal: Locked Bag 8,
Gosford NSW, Australia 2250
Téléfono: ++(61 2) 4340 2044
Fax: ++(61 2) 4340 1991
Correo electrónico: njp1@njphillips.com.au
Sitio web: www.njphillips.com

HECHO EN AUSTRALIA
Adresse Postale: Locked Bag 8,
Gosford NSW, Australia 2250
Téléphone: +612 4340 2044
Fax: +612 4340 1991
E-mail: njp1@njphillips.com.au
www.njphillips.com

PRODOTTO IN AUSTRALIA
Indirizzo postale: Locked Bag 8
Gosford NSW 2250 AUSTRALIA
TELEFONO: +61 2 4340 2044
FAX: +61 2 4340 1991
E-mail: njp1@njphillips.com.au
www.njphillips.com

QL405-R12

AS OUR POLICY IS ONE OF CONTINUOUS IMPROVEMENT, THE MANUFACTURER RESERVES THE RIGHT TO ALTER THESE SPECIFICATIONS AT ANY TIME. ALL PRODUCTS PRODUCED BY N.J. PHILLIPS PTY LIMITED, ARE IDENTIFIED BY A UNIQUE BATCH NUMBER. THIS IDENTIFICATION NUMBER IS AFFIXED TO THE PRODUCT TO ALLOW TRACEABILITY BY THE MANUFACTURER AND MUST NOT BE REMOVED IF PRODUCT INTEGRITY IS TO BE MAINTAINED.

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