

# NJ Phillips Plastic F-Grip Dial-a-Dose Injector Range

## 'F-Grip' Dial-a-Dose Bottle Attachment\*

Code: PAS1395 (1ml)

Code: PAS1397 (2ml)

Code: PAS1399 (5ml)

Quick Release  
Bottle Fitting  
for securing the  
vaccine pack.

MODEL: PAS1397

MIM Nut  
Metal Integrated  
Mount

Quick and easy  
Dial-a-Dose  
adjustment.

## 'F-Grip' Dial-a-Dose\*

Code: PAS1394 (1ml)

Code: PAS1396 (2ml)

Code: PAS1398 (5ml)


Standard 3/16"  
(4mm) inlet  
fitting, 1/4  
(6mm) option  
available.

MODEL: PAS1396

## Suitable for: Cattle, Pigs & Sheep

- For use with a wide variety of injectable solutions.
- Lightweight, easy and accurate to use.
- Superior comfort grip and easy-off needle nut.
- Supplied with 2 needles.
- Sterilisable by boiling.
- Feedtube included.

\* 10ml available on request, minimum order applies.



**njphillips.com**

*Our expertise is in your hands.*

**NJ PHILLIPS PTY LIMITED** ABN 36 000 082 002

44 Gindurra Road, Somersby NSW, Australia 2250

Postal Address: Locked Bag 8, Gosford NSW, Australia 2250

Phone: +61 2 4340 2044 Fax: +61 2 4340 1991

Email: njp1@njphillips.com.au



The information contained in this document is legally privileged and confidential intended only for the use of the individual or entity. If the receiver of this document is not the intended recipient, the receiver is hereby notified that any dissemination, distribution or copy of this information is strictly prohibited. If this information is received in error, please notify us immediately. As our policy is one of continuous improvement, the manufacturer reserves the right to alter these specifications at any time. All products produced by NJ Phillips Pty Limited are identified by a unique batch number which is affixed to the product to allow traceability by the manufacturer and must not be removed if product integrity is to be maintained. Not all part numbers may be available.