

NJ Phillips Plastic Push-Grip Injector

2ml (PAS1148B) 5ml (PAS1149B) 10ml (PAS1150B)



SXM52 Bottle Mount Attachment

PAS1149B shown with Bottle Mount fitted.

SH7 Feed Tube

PAS1149B shown without Bottle Mount fitted.

Suitable for: Cattle, Pigs & Sheep

- The bottle mount is well-balanced and easy to fit as the bottle is centred over the injector.
- 'Dial-a-dose' system provides 10 quick and easy to select doses.
- Dual-purpose injector easily changed from bottle mount to direct feed tube.
- New ergonomic *push-action* design.
- Ready for use as a bottle mount injector for standard 20mm neck vaccine bottles.
- Includes Bottle Mount Attachment (Code: SXM552). Also available separately.
- We can supply Draw-off Cap WX231 for neck sizes 28mm or 30mm.

njphillips.com

NJ PHILLIPS PTY LIMITED ABN 36 000 082 002
44 Gindurra Road, Somersby NSW, Australia 2250 Postal Address: Locked Bag 8, Gosford NSW, Australia 2250
Phone:+61 2 4340 2044 Fax:+61 2 4340 1991 Email:njp1@njphillips.com.au

Our expertise is in your hands.

The information contained in this document is legally privileged and confidential intended only for the use of the individual or entity. If the receiver of this document is not the intended recipient, the receiver is hereby notified that any dissemination, distribution or copy of this information is strictly prohibited. If this information is received in error, please notify us immediately. As our policy is one of continuous improvement, the manufacturer reserves the right to alter these specifications at any time. All products produced by NJ Phillips Pty Limited are identified by a unique batch number which is affixed to the product to allow traceability by the manufacturer and must not be removed if product integrity is to be maintained. Not all part numbers may be available.